



Recombinant Solanum tuberosum Light-induced protein, chloroplastic

Product Code	CSB-MP302156FIG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P80471
Product Type	Recombinant Protein
Immunogen Species	Solanum tuberosum (Potato)
Purity	>85% (SDS-PAGE)
Sequence	ATNYDKE DEWGPEVEQI RPPGVAVVEE EPPKEPSEIE LLKKQLADSL YGTNRGLSAS SETRAEIVEL ITQLESKNPN PAPTEAL TLL NGKWILAYTS FSGLP LLSR GNLPLVRVEE ISQTIDSESF TVQNSVVFAG PLATTSISTN AKFEVRSPKR VQIKFEEGII GTPQLTDSIV LPENVEFLGQ KIDLSPFKGL ITSVQDTASS VAKSISSQPP IKFPITNNA QSWLLTTYLD DELRISRGDA GSVFVLIKEG SPLLLKP
Source	Mammalian cell
Protein Names	Recommended name: Light-induced protein, chloroplastic Alternative name(s): Drought-induced stress protein CDSP-34
Expression Region	64-326
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.