



# Recombinant Influenza A virus Nucleoprotein (NP)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-BP302450IGF                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Uniprot No.</b>       | P69296                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Immunogen Species</b> | Influenza A virus (strain A/Shanghai/16/1989 H3N2)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Sequence</b>          | <p>MASQGTKRSY EQMETDGERQ NATEIRASVG KMIDGIGRFY IQMCTELKLS<br/> DYEGRLIQNS LTVERMVLSA FDERRNRYLE EHPSAGKDPK KTGGPIYKRV<br/> GGRWMRELVL YDKEEIRRIW RQANNGDDAT RGLTHMMIWH SNLNDTTYQR<br/> TRALVRTGMD PRMCSLMQGS TLPRRSGAAG AAVKGIGTMV MELIRMIKRG<br/> INDRNFWRGE NGRKTRSAYE RMCNILKGKF QTAAQRAMMD<br/> QVRESRNPNGN AEIEDLIFSA RSALILRGSV AHKSCLPACV YGPAVSSGYD<br/> FEKEGYSLVG IDPFKLLQNS QVYSLIRPNE NPAHKSQLVW MACHSAAFED<br/> LRLLSFIRGT KVSPRGKLST RGVQIASNEN MDNMESSTLE LRSRYWAIRT<br/> RSGGNTNQQR ASAGQISVQP TFSVQRNLPF EKSTVMAAFT GNTGRTSDM<br/> RAEIIRMMEG AKPEEVSFRRG RGVFELSDEK ATNPVPSFD MSNEGSYFFG<br/> DNAEEYDN</p> |
| <b>Source</b>            | Baculovirus                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Target Names</b>      | NP                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Protein Names</b>     | Recommended name: Nucleoprotein Alternative name(s): Nucleocapsid protein<br>Short name= Protein N                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Expression Region</b> | 1-498                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                                                                                                                                                                  |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                                                                                                                           |