



Recombinant Escherichia coli Putative colanic acid biosynthesis glycosyltransferase wcaL (wcaL)

Product Code	CSB-BP302568ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P71243
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MKVGFFLLKF PLSSETFVLN QITAFIDMGF EVEILALQKG DTQNTHAAWT KYNLAARTRW LQDEPTGKVA KLRHRASQTL RGIHRKNTWQ ALNLKRYGAE SRNLILSAIC GQVATPFRAD VFIAHFGPAG VTAACLRELG VIRGKIATIF HGIDISSREV LNHYTPEYQQ LFRRGDLMLP ISDLWAGRLQ KMGCPREKIA VSRMGVDMTR FSPRPVKAPA TPLEIISVAR LTEKKGLHVA IEACRQLKEQ GVAFRYRILG IGPWERRLRT LIEQYQLEDV VEMPGFKPSH EVKAMLDDAD VFLLPSTVGA DGDMEGIPVA LMEAMAVGIP VVSTLHSGIP ELVEADKSGW LVPENDARAL AQRLAAFSQL DTDELAPVVK RAREKVEHDF NQQVINRELA SLLQAL
Source	Baculovirus
Target Names	wcaL
Protein Names	Recommended name: Putative colanic acid biosynthesis glycosyltransferase wcaL EC= 2.4.-.-
Expression Region	1-406
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.