



Recombinant Escherichia coli Uncharacterized protein YmdC (ymdC)

Product Code	CSB-MP302795ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P75919
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MPRLASAVLP LCSQHPGQCG LFPLEKSLDA FAARYRLAEM AEHTLDVQYY IWQDDMSGRL LFSALLAAAK RGVVRVRLLD DNNTPLGDDI LRLLDHPRI EVRLFNPFSF RLLRPLGYIT DFSRLNRRMH NKSFTVDGVV TLVGGRNIGD AYFGAGEEPL FSDLVMAIG PVVEDVADDY ARYWYCKSVS PLQQVLDVPE GEMADRIELP ASWHNDAMTH RYLKRMESP FINHLVDGTL PLIWAKTRLL SDDPAKGEGK AKRHSLLPQR LFDIMGSPSE RIDISSYFV PTRAGVAQLL RMVRKGVKIA ILTNSLAAND VAVVHAGYAR WRKLLRYGV ELYELKPTRE QSSTLHDRGI TGNSGASLHA KTFSIDGKTV FIGSFNFDPR STLLNTEMGF VIESETLAQL IDKRFIQSQY DAAWQLRLDR WGRINWVDRH AKKEILKKE PATSFWKRVV VRLASILPVE WLL
Source	Mammalian cell
Target Names	clsC
Protein Names	Recommended name: Uncharacterized protein YmdC
Expression Region	1-473
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.