



Recombinant Escherichia coli Transposase insH for insertion sequence element IS5Y (insH5)

Product Code	CSB-EP302807ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P76071
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MSHQLTFADS EFSTKRRQTR KEIFLSRMEQ ILPWQNMTAV IEPFYPKAGN GRRPYPLETM LRIHCMQHWY NLSDGAMEDA LYEIASMRLF ARLSLDSALP DRTTIMNFRH LLEQHQLARQ LFKTINRWLA EAGVMMTQGT LVDATIIEAP SSTKNKEQQR DPEMHQTKKG NQWHFGMKAH IGVDAKSGLT HSLVTTAANE HDLNQLGNLL HGEEQFVSAD AGYQGAPQRE ELAEVDVDWL IAERPGKVKT LKQNPRKNKT AINIEYMKAS IRARVEHPFR IIKRQFGFVK ARYKGLLKND NQLAMLFTLA NLFRVDQMIR QWERSQ
Source	E.coli
Target Names	insH5
Protein Names	Recommended name: Transposase insH for insertion sequence element IS5Y
Expression Region	1-326
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.