



Recombinant Escherichia coli Uncharacterized protein YncE (yncE)

Product Code	CSB-MP302812ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P76116
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	≥85% (SDS-PAGE)
Sequence	AEEMLRKAVG KGAYEMAYSQ QENALWLATS QSRKLDKGGV VYRLDPVTLE VTQAIHNDLK PFGATINNTT QTLWFGNTVN SAVTAIDAKT GEVKGRVLVD DRKRTEEV RP LQPRELVADD ATNTVYISGI GKESVIWVVD GGNIKLKTAI QNTGKMSTGL ALDSEGKRLY TTNADGELIT IDTADNKILS RKKLLDDGKE HFFINISLDT ARQRAFITDS KAAEVLVVD T RNGNILAKVA APESLAVLFN PARNEAYVTH RQAGKVSVID AKSYKVVKTF DTPTHPNSLA LSADGKTL YV SVKQKSTKQQ EATQPDDVIR IAL
Source	Mammalian cell
Target Names	yncE
Protein Names	Recommended name: Uncharacterized protein YncE
Expression Region	31-353
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.