



Recombinant Escherichia coli 5-hydroxyisourate hydrolase (hiuH)

Product Code	CSB-MP302828ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P76341
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	AQQNILS VHILNQQTGK PAADVTVTLE KKADNGWLQL NTAKTDKDGR IKALWPEQTA TTGDYRVVFK TGDYFKKQNL ESFFPEIPVE FHINKVNEHY HVPILLSQYG YSTYRGS
Source	Mammalian cell
Target Names	hiuH
Protein Names	Recommended name: 5-hydroxyisourate hydrolase Short name= HIU hydrolase Short name= HIUHase EC= 3.5.2.17 Alternative name(s): Transthyretin-like protein Short name= TLP Transthyretin-related protein Short name=
Expression Region	24-137
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.