



# Recombinant Escherichia coli Uncharacterized HTH-type transcriptional regulator yeeY (yeeY)

<b>Product Code</b>	CSB-EP302830ENV-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P76369
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain K12)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKPLLDVLM I L D A L E K E G S F A A A S A K L Y K T P S A L S Y T V H K L E S D L N I Q L L D R S G H R A K F T R T G K M L L E K G R E V L H T V R E L E K Q A I K L H E G W E N E L V I G V D D T F P F S L L A P L I E A F Y Q H H S V T R L K F I N G V L G G S W D A L T Q G R A D I I V G A M H E P P S S E F G F S R L G D L E Q V F A V A P H H P L A L E E E P L N R R I I K R Y R A I V V G D T A Q A G A S T A S Q L L D E Q E A I T V F D F K T K L E L Q I S G L G C G Y L P R Y L A Q R F L D S G A L I E K K V V A Q T L F E P V W I G W N E Q T A G L A S G W W R D E I L A N S A I A G V Y A K S D D G K S A I
<b>Source</b>	E.coli
<b>Target Names</b>	yeeY
<b>Protein Names</b>	Recommended name: Uncharacterized HTH-type transcriptional regulator yeeY
<b>Expression Region</b>	1-309
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.