



Recombinant Escherichia coli Uncharacterized protein YdeU (ydeU)

Product Code	CSB-MP302857ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P77286
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	≥85% (SDS-PAGE)
Sequence	<p>MNSEGGKPGN VLTVNGNYTG NNGLMTFNAT LGGDNSPTDK MNVKGDTQGN TRVRVDNIGG VGAQTVNGIE LIEVGGNSAG NFALTTGTVE AGAYVYTLAK GKGNDKKNWY LTSKWVGTP ADTPDPINNP PVVDPEGPSV YRPEAGSYIS NIAAANSLFS HRLHDRLGEP QYTDLSHSQG SASSMWMRHV GGHERSRAGD GQLNTQANRY VLQLGGDLAQ WSSNAQDRWH LGVMAGYANQ HSNTQSNRVG YKSDGRISGY SAGLYATWYQ NDANKTGAYV DSWALYNWFD NSVSSDNRSA DDYDSRGVTA SVEGGYTFEA GTFSGSEGL NTWYVQPQAQ ITWMGVKDSD HTRKDGTRIE TEGDGNVQTR LGVKTYLNH HQRDDGKQRE FQPYIEANWI NNSKVYAVKM NGQTVGREGA RNLGEVRTGV EAKVNNNLSL WGNVGVQLGD KGYSDTQGML GVKYSW</p>
Source	Mammalian cell
Target Names	ydeU
Protein Names	Recommended name: Uncharacterized protein YdeU
Expression Region	1-466
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.