



Recombinant Escherichia coli Putative xanthine dehydrogenase YagS FAD-binding subunit (yagS)

Product Code	CSB-BP302860ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P77324
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MKAFTYERVN TPAEAALSAQ RVPGAKFIAG GTNLLDLMKL EIETPTHLID VNLGLDKIE VTDAGGLRIG ALVRNTDLAA HERVRRDYAV LSRALLAGAS GQLRNQATTA GNLLQRTRCP YFYDTNQCPCN KRLPGSGCAA LEGFSRQHAV VGVSEACIAT HPSDMAVAMR LLDVAVVETIT PEGKTRSITL ADFYHPPGKT PHIETALLPG ELIVAVTLPP PLGGKHIYRK VRDRASYAFA LVSVAIIQP DGSGRVALGG VAHKPWRIEA ADAQLSQGAQ AVYDTLFASA HPTAENTFKL LLAKRTLASV LAEARAQA
Source	Baculovirus
Target Names	yagS
Protein Names	Recommended name: Putative xanthine dehydrogenase YagS FAD-binding subunit EC= 1.17.1.4
Expression Region	1-318
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.