



Recombinant Escherichia coli Putative tail fiber assembly protein homolog from lambdoid prophage DLP12 (tfaD)

Product Code	CSB-BP302881ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P77699
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	≥85% (SDS-PAGE)
Sequence	MNISNSQVNR LRHFVFRAGLR SLFRPEPQTA VEWADANYYL PKESAYQEGR WETLPFQRAI MNAMGSDYIR EVNVVKSARV GYSKMLLGVIY AYFIEHKQRN TLIPAGFVAV FNSDESSWHL VEDHRGKTVY DVASGDALFI SELGPLPENV TWLSPEGEFQ KWNGTAWVKD AEAEKLFRIE EAEETKNSLM QVASEHIAPL QDAVDLEIAT EEETSLLLEAW KKYRVLLNRV DTSTAPDIEW PTNPVRE
Source	Baculovirus
Target Names	tfaD
Protein Names	Recommended name: Putative tail fiber assembly protein homolog from lambdoid prophage DLP12
Expression Region	1-247
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.