



Recombinant Escherichia coli Putative uncharacterized protein ybfO (ybfO)

Product Code	CSB-EP302885ENV-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P77779
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	≥85% (SDS-PAGE)
Sequence	MWPDNRIARD AHYLYRYDRH GRLTEKTDLI PEGVIRTDDE RTHRYHYDSQ HRLVHYTRTQ YAEPLVESRY LYDPLGRRVA KRVWRRERDL TGWMSLSRKP QVTWYGWDGD RLTTIQNDRT RIQTIYQPGS FTPLIRVETA TGELAKTQRR SLADTLQQSG GEDGGSVVFP PVLVQMLDRL ESEILADRVS EESRRWLASC GLTVEQMQRNQ MDPVYTPARK IHLYHCDHRG LPLALVSTEG ATEWCAEYDE WGNLLNEENP HQLQQILRLP GQQYDEESGL YYNRHRYYDP LQGRYITQDP IGLKGWVNFY QYPLNPVQYI DSMGLASKYG HLNNGGYGAR PNKPPTPDPS KLPDIAKQLR LPYPIDQASS APNVFKTFFR ALSPYDYTLY CRKWVKPNLT CTPQDDSQYP GMDTKTASDY LPQTNWPTTQ LPPGYTCAEP YLFPDINKPD GPATAGIDDL GEILAKMKQR TSRGIRK
Source	E.coli
Target Names	ybfO
Protein Names	Recommended name: Putative uncharacterized protein ybfO
Expression Region	1-477
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.