



Recombinant Enterobacteria phage 933W Protein ninG (ninG)

Product Code	CSB-EP303117ECC
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P69174
Product Type	Recombinant Protein
Immunogen Species	Enterobacteria phage 933W (Bacteriophage 933W)
Purity	>85% (SDS-PAGE)
Sequence	MAKPARRKCK ICKEWFHPAF SNQWWCCPEH GTQLALERRS KEREKAEKAA EKKRRREEQK QKDKLKIRKL ALKPRSYWIK QAQQAVNAFI RERDRDLPCI SCGTLTSAQW DAGHYRTTAA APQLRFDERN IHKQCVCNQ HKSIGNLVPYR VELINRIGQE AVDEIESNHN RHRWTVVECR AIKAKYQQKL KDLRNSRSEA A
Source	E.coli
Target Names	ninG
Protein Names	Recommended name: Protein ninG
Expression Region	1-201
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.