



Recombinant Bacteroides fragilis Elongation factor P (efp)

Product Code	CSB-BP303228BDQ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P70889
Product Type	Recombinant Protein
Immunogen Species	Bacteroides fragilis (strain YCH46)
Purity	>85% (SDS-PAGE)
Sequence	MINAQDIKNG TCIRMDGKLY FCIEFLHVKP GKGNTFMRTK LKDVVSGYVL ERRFNIGEKL EDVRVERRPY QYLYKEGEDY IFMNQETFDQ HPIAHDLING VDFLLEGAVV EVVSDASTET VLYADMPIKV QMKVITYTEPG LKGDTATNTL KPATVESGAT VRVPLFISEG ETIEIDTRDG SYVGRVKA
Source	Baculovirus
Target Names	efp
Protein Names	Recommended name: Elongation factor P Short name= EF-P
Expression Region	1-188
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.