



Recombinant Human Ubiquitin carboxyl-terminal hydrolase 46 (USP46)

| | |
|--------------------------|---|
| Product Code | CSB-EP025738HU |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P62068 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | <p>MTVRNIASICNMGTNASALEKDIGPEQFPINEHYFGLVNFNGNTCYCNSVLQALY FCRPFR ENVLAYKAQQKKKENLLTCLADLFHSIATQKKKVGVIPPKKFISRLRKENDLFDN YMQQD AHEFLNYLLNTIADILQEEKKQEKQNGKLNKNGMNEPAENNKPELTWVHEIFQ GTLTNET RCLNCETVSSKDEDFLDLSVDVEQNTSITHCLRDFSNTETLCSEQKYCETCC SKQEAQK RMRVKKLPMILALHLKRFKYMEQLHRYTKLSYRVVFPLELRLFNTSSDAVNLDLDR MYDLVA VVVHCGSGPNRGHYITIVKSHGFWLLFDDDDIVEKIDAQAIEEFYGLTSDISKNSE SGYIL FYQSRE</p> |
| Source | E.coli |
| Target Names | USP46 |
| Protein Names | <p>Recommended name: Ubiquitin carboxyl-terminal hydrolase 46 EC= 3.4.19.12 Alternative name(s): Deubiquitinating enzyme 46 Ubiquitin thioesterase 46 Ubiquitin-specific-processing protease 46</p> |
| Expression Region | 1-366 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | <p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p> |