



Recombinant Escherichia coli Uncharacterized protein ybfM (ybfM)

Product Code	CSB-MP303448ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P75733
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	EGFIDDST LTGGIYYWQR ERDRKDVTDG DKYKTNLSHS TWNANLDFQS GYAADMFGLD IAAFTAIEMA ENGDSSHPNE IAFSKSNKAY DEDWSGDKSG ISLYKAAAKF KYGPVWARAG YIQPTGQTLT APHWSFMPGT YQGAEAGANF DYGDAGALSF SYMWTNEYKA PWHLEMDEFY QNDKTTKVDY LHSFGAKYDF KNNFVLEAAF GQAEGYIDQY FAKASYKFDI AGSPLTTSYQ FYGTRDKVDD RSVNDLYDGT AWLQALTFGY RAADVVDLRL EGTWVKADGQ QGYFLQRMTP TYASSNGRLD IWWDNRSDFN ANGEKAVFFG AMYDLKNWNL PGFAIGASYV YAWDAKPATW QSNPDAYYDK NRTIEESAYS LDAVYTIQDG RAKGTMFKLH FTEYDNHSDI PSWGGGYGNI FQDERDVKFM VIAPFTIF
Source	Mammalian cell
Target Names	chiP
Protein Names	Recommended name: Uncharacterized protein ybfM
Expression Region	33-468
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.