



# Recombinant Escherichia coli Putative Lon protease homolog (ycbZ)

<b>Product Code</b>	CSB-EP303458ENV-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P75867
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain K12)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTITKLAWRD LVPDTSYQE IFAQPHLIDE NDPLFSDTQP RLQFALEQLL HTRASSSFML AKAPEESEYL NLIANAARTL QSDAGQLVGG HYEVSQHSIR LRHAVSADDN FATLTQVVAA DWVEAEQLFG CLRQFNGDIT LQPLVHQQAN GGILIISLRT LLAQPLLWMR LKNIVNRERF DWVAFDESRP LPVSVPSMPL KLVILVGER ESLADFEQEME PELSEQAIYS EFEDTLQIVD AESVTQWCRW VTFTARHNHL PAPGADAWPI LIREAARYTG EQETLPLSPQ WILRQCKEVA SLCDGDTFSG EQLNMLQQR EWREGFLAER MQDEILQEQLIETEGERIG QINALSVIEF PGHPRAFGEPRSRISCVVHIG DGEFTDIERK AELGGNIHAK GMMIMQAFLM SELQLEQQIP FSASLTFEQS YSEVDGDSAS MAELCALISA LADVPVNQSI AITGSVDQFG RAQPVGGLNE KIEGFFAICQ QRELTGKQGV IIPANVRHL SLHSELVKAV EEGKFTIWA VDDVTDALPLL LNLVWDGEGQ TTLMQTIQER IAQASQQEGR HRFWPLRWL NWFIPN
<b>Source</b>	E.coli
<b>Target Names</b>	ycbZ
<b>Protein Names</b>	Recommended name: Putative Lon protease homolog EC= 3.4.21.- Alternative name(s): ATP-dependent protease La homolog
<b>Expression Region</b>	1-586
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.