



Recombinant Escherichia coli Putative Lon protease homolog (ycbZ)

Product Code	CSB-EP303458ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P75867
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	<p>MTITKLAWRD LVPDTSYQE IFAQPHLIDE NDPLFSDTQP RLQFALEQLL HTRASSSFML AKAPEESEYL NLIANAARTL QSDAGQLVGG HYEVSQHSIR LRHAVSADDN FATLTQVVAA DWVEAEQLFG CLRQFNGDIT LQPLVHQQAN GGILIISLRT LLAQPLLWMR LKNIVNRERF DWVAFDESRP LPVSVPSMPL KLVILVGER ESLADFQEME PELSEQAIYS EFEDTLQIVD AESVTQWCRW VTFTARHNHL PAPGADAWPI LIREAARYTG EQETLPLSPQ WILRQCKEVA SLCDGDTFSG EQLNMLQQR EWREGFLAER MQDEILQEQLIETEGERIG QINALSVIEF PGHPRAFGEPR SRISCVVHIG DGEFTDIERK AELGGNIHAK GMMIMQAFLM SELQLEQQIP FSASLTFEQS YSEVDGDSAS MAELCALISA LADVPVNQSI AITGSVDQFG RAQPVGGLNE KIEGFFAICQ QRELTGKQGV IIPATANVRHL SLHSELVKAV EEGKFTIWAV DDVTDALPLL LNLVWDGEGQ TTLMQTIQER IAQASQQEGR HRFWPLRWL NWFIPN</p>
Source	E.coli
Target Names	ycbZ
Protein Names	Recommended name: Putative Lon protease homolog EC= 3.4.21.- Alternative name(s): ATP-dependent protease La homolog
Expression Region	1-586
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.