



Recombinant Escherichia coli Putative ABC transporter periplasmic-binding protein ycjN (ycjN)

Product Code	CSB-EP303471ENV-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P76042
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	G CKEENKTNVS IEFMHSSVEQ ERQAVISKLI ARFEKENPGI TVKQVPVEED AYNTKVITLS RSGSLPEVIE TSHDYAKVMD KEQLIDRKAV ATVISNVGEG AFYDGVLRIV RTEDGSAWTG VPVSAWIGGI WYRKDVLAKE GLEPKNWQQ LLDVAQKLND PANKKYGIAL PTAESVLTEQ SFSQFALSQ ANVFNAEGKI TLDTPEMMQA LTYYRDLTAN TMPGSNDIME VKDAFMNGTA PMAIYSTYIL PAVIKEGDPK NVGFVVPTEK NSAVYGMLTS LTITAGQKTE ETEAAEFVVT FMEQADNIAD WMMSPGAAL PVNKAVTTA TWKDNDVIKA LGELPNQLIG ELPNIQVFGA VGDKNFTRMG DVTGSGVVSS MVHNVTVGKA DLSTTLQASQ KKLDELIEQH
Source	E.coli
Target Names	ycjN
Protein Names	Recommended name: Putative ABC transporter periplasmic-binding protein ycjN
Expression Region	20-430
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.