



Recombinant Escherichia coli Uncharacterized electron transport protein ykgF (ykgF)

Product Code	CSB-BP303539ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P77536
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MSIKTSNTDF KTRIRQQIED PIMRKAVANA QQRIGANRQK MVDELGHWEE WRDRAAQIRD HVLSNLDAYL YQLSEKVTQN GGHVYFARTK EDATRYILQV AQRKNARKVV KSKSMVTEEI GVNHLVLDAG IQVIETDLGE YILQLDQDPP SHVVVPAIHK DRHQIRRVLH ERLGYEGPET PEAMTLFIRQ KIREDFLSAE IGITGCNFAV AETGSVCLVT NEGNARMCTT LPKTHIAVMG MERIAPTFAE VDVLITMLAR SAVGARLTGY NTWLTGPREA GHVDGPEEFH LVIVDNGRSE VLASEFRDVL RCIRCGACMN TCPAYRHIGG HGYGSIYPGP IGAVISPLLG GYKDFKDLPY ACSLCTACDN VCPVRIPLSK LILRHRVMA EKGITAKAEQ RAIKMFAYAN SHPGLWKVGM MAGAHAASWF INGGKTPLKF GAISDWMEAR DLPEADGESF RSWFKKHQAQ EKKNQ
Source	Baculovirus
Target Names	ykgF
Protein Names	Recommended name: Uncharacterized electron transport protein ykgF
Expression Region	1-475
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.