



Recombinant Escherichia coli Long-chain acyl-CoA thioesterase FadM (fadM)

Product Code	CSB-YP303548ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P77712
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	≥85% (SDS-PAGE)
Sequence	MQTQIKVIRGY HLDVYQHVNN ARYLEFLEEA RWDGLENLSDS FQWMTAHNIA FVVVNININY RRPVAVLSDLL TITSQLQQLN GKSGILSQVI TLEPEGQVVA DALITFVCID LKTQKALALE GELREKLEQM VK
Source	Yeast
Target Names	fadM
Protein Names	Recommended name: Long-chain acyl-CoA thioesterase FadM EC= 3.1.2.- Alternative name(s): Acyl-CoA thioester hydrolase Thioesterase 3
Expression Region	1-132
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.