



Recombinant 3-deoxy-D-manno-octulosonate 8-phosphate phosphatas KdsC (kdsC)

Product Code	CSB-EP303673EGX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P67653
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MSKAGASLAT CYGPVSADVM AKAENIRLLI LDVDGVLSDG LIYMGNNGEE LKAFNVRDGY GIRCALTS DI EVAITGRKA KLVEDRCATL GITHLYQQQS NKLIAFSDLL EKLAIAPENV AYVGDDLIDW PVMEKVGLSV AVADAHPLLI PRADYVTRIA GGARGAVREVC DLLLLAQQKL DEAKGQSI
Source	E.coli
Target Names	kdsC
Protein Names	Recommended name: 3-deoxy-D-manno-octulosonate 8-phosphate phosphatas KdsC EC= 3.1.3.45 Alternative name(s): KDO 8-P phosphatase
Expression Region	1-188
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.