



Recombinant Rice stripe virus Nucleoprotein (pc3)

Product Code	CSB-BP303737REQ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P68559
Product Type	Recombinant Protein
Immunogen Species	Rice stripe virus (isolate T) (RSV)
Purity	>85% (SDS-PAGE)
Sequence	MGTNKPATLA DLQKAINDIS KDALSYLTAH KADVVTFAGQ IEYAGYDAAT LIGILKDKGG DTLAKDMTMC ITMRYVRGTG FVRDVTKKVK VAAGSTEAST LVSRYGIVSS VGTNANAILL GRLAQLFPNV SHEVVRQISG VKMAVDSSDL GLTGCDNLLW DYVPQYIKLE SETAPYCSTH SLSHILFVVH IHSFQITKK TMPEGKKKER GLTKDIDMMK YTTGLLVITC KSKNLSDKKK EGRKKVLDE FITNGKVKTT IFDALAGMSV NTISTYGNQT RLYLAQQSKL MKILAENTSK TATEVSGLVK EFFEDEAEGA DD
Source	Baculovirus
Target Names	pc3
Protein Names	Recommended name: Nucleoprotein Alternative name(s): Coat protein Short name= CP Nucleocapsid protein Short name= Protein N Protein pc3
Expression Region	1-322
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.