



Recombinant Human Eukaryotic translation initiation factor 3 subunit E (EIF3E)

Product Code	CSB-MP007534HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P60228
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAEYDLTTRI AHFLDRHLVF PLLEFLSVKE IYNEKELLQG KLDLLSDTNM VDFAMDVYKN LYSDDIPHAL REKRTTVVAQ LKQLQAETEP IVKMFEDPET TRQMQSTRDG RMLFDYLADK HGFRQEYLDL LYRYAKFQYE CGNYSGAAEY LYFFRVLVPA TDRNALSSLW GKLASEILMQ NWDAAMEDLT RLKETIDNNS VSSPLQSLQQ RTWLIHWSLF VFFNHPKGRD NIIDLFLYQP QYLNAIQTMC PHILRYLTTA VITNKDVRKR RQVLKDLVKV IQQESYTYKD PITEFVECLY VNFDFDGAQK KLRECESVLV NDFFLVACLE DFIENARLFI FETFCRIHQ ISINMLADKL NMTPEEAERW IVNLIRNARL DAKIDSKLGH VVMGNNAVSP YQQVIEKTKS LSFERSQMLAM NIEKLNQNS RSEAPNWATQ DSGFY
Source	Mammalian cell
Target Names	EIF3E
Protein Names	Recommended name: Eukaryotic translation initiation factor 3 subunit E Short name= eIF3e Alternative name(s): Eukaryotic translation initiation factor 3 subunit 6 Viral integration site protein INT-6 homolog eIF-3 p48
Expression Region	1-445
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.