



Recombinant Escherichia coli Iron-sulfur cluster repair protein YtfE (ytfE)

Product Code	CSB-BP303998ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P69506
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MAYRDQPLGE LALSIPRASA LFRKYDMDYC CGGKQTLARA AARKELDVEV IEAELAKLAE QPIEKDWRSA PLAEIIDHII VRYHDRHREQ LPELILQATK VERVHADKPS VPKGLTKYLT MLHEELSSHM MKEEQILFPM IKQGMGMSQAM GPISVMESEH DEAGELLEVI KHTTNNVTPP PEACTIONWKAM YNGINELIDD LMDHISLENN VLFPRALAGE
Source	Baculovirus
Target Names	ytfE
Protein Names	Recommended name: Iron-sulfur cluster repair protein YtfE Alternative name(s): Regulator of cell morphogenesis and NO signaling Short name= RCMNS
Expression Region	1-220
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.