



Recombinant *Acidianus ambivalens* Signal recognition 54 kDa protein (srp54)

Product Code	CSB-MP304063ABO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P70722
Product Type	Recombinant Protein
Immunogen Species	<i>Acidianus ambivalens</i> (<i>Desulfurolobus ambivalens</i>)
Purity	≥85% (SDS-PAGE)
Sequence	SKLLDNLRDA VRKFLTGSSS YDKAVEDFIK ELQKSLISAD VNVKLVFSLT NKIKERLKNE KPPTYIERRE WFIKIVYDEL SNLFGGDKEP KVIPDKIPYV IMLVGVQGTG KTTTAGKLAY FYKKKGFKVG LVGADVYRPA ALEQLQQLGQ QIGVPVYGEP GEKDAVGIK RGVEKFLSEK MEIIIVDTAG RHGYGEEAAL LEEMKNIYEA IKPDEVTLVI DASIGQKAYD LASKFNQASK IGTIITKMD GTAKGGGALS AVAATGATIK FIGTGEKIDE LEVFNPRRFV ARILGMGDIE TILEKIKEVE NYDKMQKKME EVISGKGKLT LRDVYNQLIA LRKMGPLSKL FQLLPGIGML GQIPEDQLKV GEEKMRKWLA IMNSMTYEEL DNPSIIDKSR MRRIALGSGT EIEDVKELIE HYNLMQRTIK MLKRRKKDVE KLLGQFGGES T
Source	Mammalian cell
Target Names	srp54
Protein Names	Recommended name: Signal recognition 54 kDa protein Alternative name(s): SRP54
Expression Region	1-451
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.