



# Recombinant Escherichia coli Citrate lyase alpha chain (citF)

<b>Product Code</b>	CSB-EP304297ENV-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P75726
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain K12)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTQKIEQSQR QERVAAWNRR AECDLAAFQN SPKQTYQAEK ARDRKLCANL EEAIRRSLGQ DGMTVSFHHA FRGGDLTVNM VMDVIAKMGF KNLTASSSL SDCHAPLVEH IRQGVVTRIIY TSGLRGPLAE EISRGLLAEP VQIHSHGGRV HLVQSGELNI DVAFLGVPSC DEFGNANGYT GKACCGSLGY AIVDADNAKQ VVMLTEELLP YPHNPASIEQ DQVDLIVKVD RVGDAAKIGA GATRMITTNP ELLIARSAAD VIVNSGYFKE GFSMQTGTGG ASLAVTRFLE DKMRSRDIRA DFALGGITAT MVDLHEKGLI RKLDDVQSFD SHAAQSLARN PNHIEISANQ YANWGSKGAS VDRLDVVVLS ALEIDTQFNV NVLTGSDGVL RGASGGHCDT AIASALSIIV APLVRGRIPT LVDNVLTCIT PGSSVDILVT DHGIAVNPAR PELAERLQEA GIKVVSIEWL RERARLLTGE PQPIEFTDRV VAVVRYRDGS VIDVVHQVKE
<b>Source</b>	E.coli
<b>Target Names</b>	citF
<b>Protein Names</b>	Recommended name: Citrate lyase alpha chain Short name= Citrase alpha chain EC= 4.1.3.6 Alternative name(s): Citrate (pro-3S)-lyase alpha chain Citrate CoA-transferase subunit EC= 2.8.3.10
<b>Expression Region</b>	1-510
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.