



Recombinant Escherichia coli Putative protein ybfL (ybfL)

Product Code	CSB-BP304298ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P75741
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MELKKLMEHI SIIPDYRQAW KVEHKLSGIL LLTIFAVISG AESWEDIEDF GETHLDFLKQ YGDFENGIPV HDTIARVVSC ISPAKFHECF INWMRDCHSS NDKDVAIDG KTLRHSYDKS RRRGAIHVIS AFSTMHSLVI GQIKTDKKS N EITAIPELLN MLDIKGKIIT TDAMGCQKDI AEKIQKQGGD YLFAVKGNQG RLNKAFEEKF PLKELNNPEH DSYAMSEKSH GREEIRLHIV CDVPDELIDF TFEWKGLKKL CVAVSFRSII AEQKKEPEMT VRYYYISSADL TAEKFATAIR NHWHVENKLH WRLDVVMNED DCKIRRGNAE ELFSGIRHIA INILTNDKVF KAGLRKMRK AAMDRNYLAS VLAGSGLS
Source	Baculovirus
Target Names	ybfL
Protein Names	Recommended name: Putative protein ybfL Alternative name(s): H repeat-associated protein in rhsC-phrB intergenic region ORF-H2
Expression Region	1-378
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.