



Recombinant Escherichia coli Low specificity L-threonine aldolase (ItaE)

Product Code	CSB-EP304304ENV-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P75823
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MIDLRSDTVT RPSRAMLEAM MAAPVGDDVY GDDPTVNALQ DYAAELSGKE AAIFLPTGTQ ANLVALLSHC ERGEEYIVGQ AAHNYLFEAG GAAVLGSIQP QPIDAAADGT LPLDKVAMKI KPDDIHFART KLLSLENTHN GKVLPREYLK EAWEFTRERN LALHVDGARI FNAVVAYGCE LKEITQYCD S FTICLSKGLG TPVGSLLVGN RDYIKRAIRW RKMTGGGMRQ SGILAAAGIY ALKNNVARLQ EDHDNAAWMA EQLREAGADV MRQDTNMLFV RVGEENAAAL GEYMKARNVL INASPIVRLV THLDVSREQL AEVAHWRAF LAR
Source	E.coli
Target Names	ItaE
Protein Names	Recommended name: Low specificity L-threonine aldolase Short name= Low specificity L-TA EC= 4.1.2.48
Expression Region	1-333
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.