



Recombinant Escherichia coli Gamma-glutamyl-gamma-aminobutyrate hydrolase PuuD (puuD)

Product Code	CSB-EP304321ENV-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P76038
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MENIMNPNVI GVVMCRNRLK GHATQTLQEK YLNAIHAGG LPIALPHALA EPSLLEQLLP KLDGIYLP GS PSNVQPHLYG ENGDEPDADP GRDLLSMAL NAALERRIPI FAICRGLQEL VVATGGSLHR KLCEQPELLE HREDPELPVE QQYAPSHEVQ VEEGGLLSAL LPECSNFWVN SLHGQGAKVV SPRLRVEARS PDGLVEAVSV INHPFALGVQ WHPEWNSSEY ALSRILFEGF ITACQHIAE KQRL
Source	E.coli
Target Names	puuD
Protein Names	Recommended name: Gamma-glutamyl-gamma-aminobutyrate hydrolase PuuD Short name= Gamma-Glu-GABA hydrolase EC= 3.5.1.94
Expression Region	1-254
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.