



Recombinant Escherichia coli Uncharacterized protein ydjO (ydjO)

Product Code	CSB-YP304337ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P76210
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MFSTNLIQSN YGDLNIKSLA FDSFKERLQS TMTALTFFIS TGQCDCDEIA ESNFNYMIAY MSNINYDASK PGAPALSFDY YLQDNV KYRV IINNLYGSEI RIRGINKDFI GMDVTSVFRP EKMTNLISIK NRLVIHYLRT IYEQNYIHP VGSIFAAIQK NESLLKFPSI ISMLNINLLF NPLNLP GMGS GVLEDIMSIP DSSLRKRLGY EVLSFSLQAH SLSQECIDKL DIFFADDFK YESVCIAAME HLKSKATAPI QNGPLPA
Source	Yeast
Target Names	ydjO
Protein Names	Recommended name: Uncharacterized protein ydjO
Expression Region	1-267
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.