



Recombinant Bovine Carbonic anhydrase 2 (CA2)

Product Code	CSB-MP004370BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P00921
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	>85% (SDS-PAGE)
Sequence	SHHWGYGKH NGPEHWHKDF PIANGERQSP VDIDTKAVVQ DPALKPLALV YGEATSRMV NNGHSFNVEY DDSQDKAVLK DGPLTGTYRL VQFHFHWGSS DDQGSEHTVD RKKYAAELHL VHWNTKYGDF GTAAQQPDGL AVVGVFLKVG DANPALQKVL DALDSIKTKG KSTDFPNFDP GSELLPNVDY WTYPGSLTTP PLESVTWIV LKEPISVSSQ QMLKFRTLNF NAEGPELLM LANWRPAQPL KNRQVRGFPK
Source	Mammalian cell
Target Names	CA2
Protein Names	Recommended name: Carbonic anhydrase 2 EC= 4.2.1.1 Alternative name(s): Carbonate dehydratase II Carbonic anhydrase II Short name= CA-II
Expression Region	2-260
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	CA2 is one of several (at least 7) isozymes of carbonic anhydrase. Carbonic anhydrase catalyzes reversible hydration of carbon dioxide. Defects in this enzyme are associated with osteopetrosis and renal tubular acidosis.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.