



Recombinant Escherichia coli Uncharacterized HTH-type transcriptional regulator ybdO (ybdO)

Product Code	CSB-BP304396ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P77746
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MANLYDLKKF DLNLLVIFEC IYQHLSISKA AESLYITPSA VSQSLQRLRA QFNDPLFIRS GKGIAPTTTG LNLHHHLEKN LRGLEQTINI VNKSELKKNF IYGPQLISC SNNSMLIRCL RQDSSVEIEC HDILMSAENA EELLVHRKAD LVITQMPVIS RSVICMPLHT IRNTLICSNR HPRITDNSTY EQIMAEFTQ LISKSAGVDD IQMEIDERFM NRKISFRGSS LLTIINSIAV TDLLGIVPYE LYNSYRDFLN LKEIKLEHPL PSIKLYISYN KSSLNNLVFS RFIDRLNESF
Source	Baculovirus
Target Names	ybdO
Protein Names	Recommended name: Uncharacterized HTH-type transcriptional regulator ybdO
Expression Region	1-300
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.