



Recombinant Neosartorya fumigata Major allergen Asp f 2 (AFUA_4G09580)

Product Code	CSB-MP304434NGS
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P79017
Product Type	Recombinant Protein
Immunogen Species	Neosartorya fumigata (strain ATCC MYA-4609 / Af293 / CBS 101355 / FGSC A1100) (Aspergillus fumigatus)
Purity	>85% (SDS-PAGE)
Sequence	<p> LVAT LPTSPVPIAA RATPHEPVFF SWDAGAVTSF PIHSSCNATQ RRQIEAGLNE AVELARHAKA HILRWGNESE IYRKYFGNRP TMEAVGAYDV IVNGDKANVL FRCDNPDGNC ALEGWGGHWR GANATSETVI CDRSYTTRRW LVSMCSQGYT VAGSETNTFW ASDLMHRLYH VPAVGQGWVD HFADGYDEVI ALAKSNGTES THDSEALQYF ALEAYAFDIA APGVGCAGES HGPDQGHDTG SASAPASTST SSSSSGSGSG ATTTPTDSPS ATIDVPSNCH THEGGQLHCT </p>
Source	Mammalian cell
Target Names	AFUA_4G09580
Protein Names	Recommended name: Major allergen Asp f 2 Alternative name(s): Allergen Asp f II Allergen= Asp f 2
Expression Region	17-310
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.