



# Recombinant Methanothermobacter marburgensis Putative thymidylate synthase (thyA)

<b>Product Code</b>	CSB-EP304474MSQ
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P80305
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Methanothermobacter marburgensis (strain ATCC BAA-927 / DSM 2133 / JCM 14651 / NBRC 100331 / OCM 82 / Marburg) (Methanobacterium thermoautotrophicum)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>           AYEISVDEI AEGWLKLVK IMSDGREIRD ERGLSTREVM NTVVTIKNPL            GRSGDFYHLP ARSLINIRVP EGYFWSGEKL EKYSEQFLSD DRKGFVYTYG            NRLRAHFGVD QVDRAIERLK NCKESRRATM VTWDPKIDTE SDEVPCMILV            DFKVREGRLF TTALWRSHDI YGAWFPNAV G LAYLADHVAS EVGVEVGHIT            IHSISAHIE VNFKEAKEVI KNG         </p>
<b>Source</b>	E.coli
<b>Target Names</b>	thyA
<b>Protein Names</b>	Recommended name: Putative thymidylate synthase Short name= TS Short name= TSase EC= 2.1.1.-
<b>Expression Region</b>	2-223
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.