



Recombinant Lepidoglyphus destructor Mite group 2 allergen Lep d 2

Product Code	CSB-BP304483LEA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P80384
Product Type	Recombinant Protein
Immunogen Species	Lepidoglyphus destructor (Storage mite) (Glycyphagus destructor)
Purity	>85% (SDS-PAGE)
Sequence	GKMT FKDCGHGEVT ELDITGCSGD TCVIHRGEKM TLEAKFAANQ DTAKVTIKVL AKVAGTTIQV PGLTDGCKF IKCPVKKGEA LDFIYSGTIP AITPKVKADV TAEIGDHGV MACGTVHGQV E
Source	Baculovirus
Protein Names	Recommended name: Mite group 2 allergen Lep d 2 Alternative name(s): Allergen Lep d 1 Allergen Lep d I Allergen= Lep d 2
Expression Region	17-141
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.