



Recombinant *Drosophila sechellia* Regulatory protein zeste (z)

Product Code	CSB-MP304678DMG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P67855
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila sechellia</i> (Fruit fly)
Purity	≥85% (SDS-PAGE)
Sequence	TAEKEVLYT LFHLHEEVID IKHRKKQRNK YSVRETWDKI VKDFNSHPHV SAMRNIKQIQ KFWLNSRLRK QYPYRDGSSS NLSSGSAKIS SVSVSVASAV PQQQQQQHHQ QHDNVKEEPE YQISPDASEH NPQADTFDEI EMDANDVSEI DEDPMEQQQQ QQQEAQAQAQ AQAQVQSAAA EMQKMQQVNA VAVAAAAAAN ATMINTHQIN VQISAEKLT LNDLLHFKTA RPREEIILQI KHPSEATATQ IHTIPTQAQQ HPMATITA
Source	Mammalian cell
Target Names	z
Protein Names	Recommended name: Regulatory protein zeste
Expression Region	1-268
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.