



# Recombinant Ribonucleoside-diphosphate reductase 1 subunit beta (nrdB)

<b>Product Code</b>	CSB-EP304877EOD
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P69925
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O157:H7
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>AYTTFSQTK NDQLKEPMFF GQPVNARYD QQKYDIFEKL IEKQLSFFWR  PEEVDVSRDR IDYQALPEHE KHIFISNLKY QTLDSIQGR SPNVALLPLI  SIPELETWVE TWAFSETIHS RSYTHIIRNI VNDPSVVFDD IVTNEIQKR  AEGISSYYDE LIEMTSYWHL LGEGTHTVNG KTVT VSLREL KKKLYLCLMS  VNALEAIRFY VSFACSFAPA ERELMEGNAK IIRLIARDEA LHGTGTQHML  NLLRSGADDP EMAEIAEACK QECYDLFVQA AQQEKDWADY LFRDGSMIGL  NKDILCQYVE YITNIRMQAV GLDLPFQTRS NPIPWINTWL VSDNVQVAPQ  EVEVSSYLVG QIDSEVDTDD LSNFQL</p>
<b>Source</b>	E.coli
<b>Target Names</b>	nrdB
<b>Protein Names</b>	Recommended name: Ribonucleoside-diphosphate reductase 1 subunit beta EC= 1.17.4.1 Alternative name(s): Protein B2 Protein R2 Ribonucleotide reductase 1
<b>Expression Region</b>	2-376
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.