



Recombinant Human Hemoglobin subunit delta (HBD)

Product Code	CSB-MP010152HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P02042
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	VHLTPEEKT AVNALWGKVN VDAVGGEALG RLLVVYPWTQ RFFESFGDLS SPDAVMGNPK VKAHGKVLG AFSDGLAHLN NLKGTFSQLS ELHCDKLHVD PENFRL LGNV LVCVLARNFG KEFTPQMCAA YQKVVAGVAN ALAHKYH
Source	Mammalian cell
Target Names	HBD
Protein Names	Recommended name: Hemoglobin subunit delta Alternative name(s): Delta-globin Hemoglobin delta chain
Expression Region	2-147
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	The delta (HBD) and beta (HBB) genes are normally expressed in the adult: two alpha chains plus two beta chains constitute HbA, which in normal adult life comprises about 97% of the total hemoglobin. Two alpha chains plus two delta chains constitute HbA-2, which with HbF comprises the remaining 3% of adult hemoglobin. Five beta-like globin genes are found within a 45 kb cluster on chromosome 11 in the following order: 5 -epsilon--Ggamma--Agamma--delta--beta-3 . Mutations in the delta-globin gene are associated with beta-thalassemia.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.