



Recombinant Escherichia coli Glutathione import ATP-binding protein GsiA (gsiA)

Product Code	CSB-EP305122ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P75796
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	≥85% (SDS-PAGE)
Sequence	<p>MPHSDELDAAG NVLAVENLNI AFMQDQQKIA AVRNLFSFLQ RGETLAIVGE SGSGKSVTAL ALMRLEQAG GLVQCCKMLL QRRSREVIEL SEQNAAQMRH VRGADMAMIF QEPMTSLNPV FTVGEQIAES IRLHQNASRE EAMVEAKRML DQVRIPEAQT ILSRYPHQLS GGMRQRMIA MALSCRPAVL IADEPTTALD VTIQAQILQL IKVLQKEMSM GVIFITHDMG VVAEIADRVL VMYQGEAVET GTVEQIFHAP QHPYTRALLA AVPQLGAMKG LDYPRRFPLI SLEHPAKQAP PIEQKTVDVG EPVLRVRNLV TRFPLRSGLL NRVTVREHAV EKVSFDLWPG ETLSLVGESG SGKSTTGRAL LRLVESQGGG IIFNGQRIDT LSPGKLQALR RDIQFIFQDP YASLDPQRQTI GDSIIPLRV HGLLPKGDAA ARVAWLLERV GLLPEHAWRY PHEFSGGQRQ RICIARALAL NPKVIIADEA VSALDVSIRG QIINLLDLQ RDFGIAYLFI SHDMAVVERI SHRVAVMYLG QIVEIGPRRA VFENPQHPYT RKLLAAVPVA EPSRQRQPRV LLSDDLPSNI HLRGEEVAAV SLQCVGPGHY VAQPQSEYAF MRR</p>
Source	E.coli
Target Names	gsiA
Protein Names	Recommended name: Glutathione import ATP-binding protein GsiA EC= 3.6.3.-
Expression Region	1-623
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.