



Recombinant Escherichia coli Ribosomal RNA large subunit methyltransferase I (rlmI)

| | |
|--------------------------|---|
| Product Code | CSB-EP305129ENV |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P75876 |
| Product Type | Recombinant Protein |
| Immunogen Species | Escherichia coli (strain K12) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MSVRLVLAKG REKSLRRHP WVFSGAVARM EGKASLGETI DIVDHQGKWL ARGAYSPASQ IRARVWTFDP SESIDIAFFS RRLQQAQKWR DWLAQKDGLD SYRLIAGESD GLPGITIDRF GNFLVLQLLS AGAEYQRAAL ISALQTLYPE CSIYDRSDVA VRKKEGMELT QGPVTGELPP ALLPIEEHGM KLLVDIQHGH KTGYLDQRD SRLATRYVE NKRVLNCFYSY TGGFAVSALM GGCSQVVSVD TSQEALDIAR QNVELNKLDL SKAEFVRDDV FKLLRTRYRDR GEKFDVIVMD PPKFVENKSQ LMGACRGYKD INMLAIQLLN EGGILLTFSC SGLMTSDLFQ KIIADAAIDA GRDVQFIEQF RQAADHPVIA TYPEGLYLKG FACRVM |
| Source | E.coli |
| Target Names | rlmI |
| Protein Names | Recommended name: Ribosomal RNA large subunit methyltransferase I EC=2.1.1.191 Alternative name(s): 23S rRNA m5C1962 methyltransferase rRNA (cytosine-C(5)-)-methyltransferase RlmI |
| Expression Region | 1-396 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |