



Recombinant Escherichia coli Flagellar hook protein flgE (flgE)

Product Code	CSB-MP305132ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P75937
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	≥85% (SDS-PAGE)
Sequence	AFSQAVSGL NAAATNLDVI GNNIANSATY GFKSGTASFA DMFAGSKVGL GVKVAGITQD FTDGTTTNTG RGLDVAISQN GFFRLVDSNG SVFYSRNGQF KLDENRNLVN MQGLQLTGYP ATGTPPTIQQ GANPTNISIP NTLMAAKTTT TASMQINLNS SDPLPTVTPF SASNADSYNK KGSVTVFDSQ GNAHDMSVYF VKTGDNNWQV YTQDSSDPNS IAKTATTLEF NANGTLVDGA MANNIATGAI NGAEPATFSL SFLNSMQQNT GANNIVATTQ NGYKPGDLVS YQINDDGTVV GNYSNEQTQL LGQIVLANFA NNEGLASEGD NVWSATQSSG VALLGTAGTG NFGTLTNGAL EASNVDSKE LVNMIVAQRN YQSNAQTIKT QDQILNLTVN LR
Source	Mammalian cell
Target Names	flgE
Protein Names	Recommended name: Flagellar hook protein flgE
Expression Region	2-402
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.