



Recombinant Escherichia coli Flagellar brake protein YcgR (ycgR)

Product Code	CSB-BP305139ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P76010
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MSHYHEQFLK QNPLAVLGVL RDLHKAAIPL RLSWNGGQLI SKLLAITPDK LVLDFGSQAE DNI AVLKAQH ITITAETQGA KVEFTVEQLQ QSEYLQLPAF ITVPPPTLWF VQRRRYFRIS APLHPPYFCQ TKLADNSTLR FRLYDLSLGG MGALLETAKP AELQEGMRFA QIEVNMGQWG VFHFDAQ LIS ISERKVIDGK NETITTPRLS FRFLNVSPTV ERQLQRIIFS LEREAREKAD KVRD
Source	Baculovirus
Target Names	ycgR
Protein Names	Recommended name: Flagellar brake protein YcgR Alternative name(s): Cyclic di-GMP binding protein YcgR
Expression Region	1-244
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.