



Recombinant Escherichia coli PTS-dependent dihydroxyacetone kinase operon regulatory protein (dhaR)

| | |
|--------------------------|---|
| Product Code | CSB-EP305140ENV-B |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P76016 |
| Product Type | Recombinant Protein |
| Immunogen Species | Escherichia coli (strain K12) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MSGAFNNDGR GISPLIATSW ERCNKLMKRE TWNVPHQAQG VTFASIYRRK KAMLTGQAA LEDAWYMAP RECALFILDE TACILSRNGD PQLTQQLSAL GFNDGTYCAE GIIGTCALSL AAISGQAVKT MADQHFKQVL WNWAFCATPL FDSKGRLTGT IALACPVEQT TAADLPLTLA IAREVGNLLL TDSLLAETNR HLNQLNALLE SMDDGVISWD EQGNLQFINA QAARVLRLDA TASQGRAITE LLTLPVLQQ AIKQAHPLKH VEATFESQHQ FIDAVITLKP IIETQGTSTFI LLLHPVEQMR QLMTSQLGKV SHTFAHMPQD DPQTRRLHF GRQAARSSFP VLLCGEEGVG KALLSQAIHN ESERAAGPYI AVNCELYGDA ALAEFIGGD RTDNENGRSL RLELAHGGTL FLEKIEYLAV ELQSALLQVI KQGVITRLDA RRLIPIDVKV IATTTADLAM LVEQNRFSRQ LYYALHAFEI TIPPLRMRRG SIPALVNNKL RSLEKRFSTR LKIDDDALAR LVSCAWPGND FELYSVIENL ALSSDNGRIR VSDLPEHLFT EQATDDVSAT RLSTSLFAE VEKEAIINAA QVTGGRIQEM SALLGIGRTT LWRKMKQHGI DAGQFKRRV |
| Source | E.coli |
| Target Names | dhaR |
| Protein Names | Recommended name: PTS-dependent dihydroxyacetone kinase operon regulatory protein |
| Expression Region | 1-639 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. |



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.