



Recombinant Escherichia coli Uncharacterized protein yciW (yciW)

Product Code	CSB-MP305141ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P76035
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MEQRHITGKS HWYHETQSST TEYDVLPLVP EAAKVSDPFL LDVILEKETL APFLSWLDPA RVLAVDLFPD QLTVTRSQTF TAYERLSTAL TVAQVCGVQR LCNYYSARLT PLPGPDSTRE SNHRLAQITQ YARQLASSPS IIDNRSRQHL NDVGLTAWDC VIISQIIGFI GFQARTIATF QAYLGHPVRW LPGLEIQNYA DASLFADESL RWRSSYEVEK LPEEHTKSST AELCQLAEIL SLHPISLSLL EKLLNSTRGN TQPDNQLAAL LCARINGSIPA CFATCMDSSN EYKKISTLMR KGENEINQWA DRHSVERATV QAIQWLTRAP DRFSAAQFSP LLEHEKSSTQ IINLLVWSGL CGWINRLKIA LGETY
Source	Mammalian cell
Target Names	yciW
Protein Names	Recommended name: Uncharacterized protein yciW
Expression Region	1-375
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.