



Recombinant Escherichia coli Protein ynjB (ynjB)

Product Code	CSB-BP305157ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P76223
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MRHCGWLLGL LSLFSLATHA SDWQEIKNEA KGQTVWFNAW GGDTAINRYL DWVSGEMKTH YAINLKIVRL ADAADAVKRI QTEAAAGRKT GGSVDLLWVN GENFRTLKEA NLLQTGWAET LPNWRVYDTQ LPVREDFSV TQGAESPWGG AQLTFIARRD VTPQPPQTPQ ALLEFAKANP GTVTYPRPPD FTGTAFLEQL LIMLTPDPAA LKEAPDDATF ARVTAPLWQY LDVLHPYLWR EGKDFPPSPA RMDALLKAGT LRLSLTFNPA HAQQKIASGD LPASSYSFGF REGMIGNVHF VTIPANANAS AAKVVANFL LSPDAQLRKA DPAVWGDPSV LDPQKLPDGQ RESLQSRMPQ DLPPVLAEPH AGWVNALEQE WLHRYGTH
Source	Baculovirus
Target Names	ynjB
Protein Names	Recommended name: Protein ynjB
Expression Region	1-388
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.