



# Recombinant Escherichia coli Multidrug resistance protein MdtA (mdtA)

<b>Product Code</b>	CSB-YP305167ENV
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P76397
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain K12)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	AFWFWQGRN DSRSAAPGAT KQAQQSPAGG RRGMRSGPLA PVQAATAVEQ AVPRYLTLGLG TITAANTVTV RSRVDGQLIA LHFQEGQQVK AGDLLAEIDP SQFKVALAQA QGQLAKDKAT LANARRDLAR YQQLAKTNLV SRQELDAQQA LVSETEGTIK ADEASVASAQ LQLDWSRITA PVDGRVGLKQ VDVGNQISSG DTTGIVVITQ THPIDLVFTL PESDIATVVQ AQKAGKPLVV EAWDRTNSKK LSEGTLISLD NQIDATTGTI KVKARFNNQD DALFPNQFVN ARMLVDTEQN AVVIPTAALQ MGNEGHFVWV LNSENKVSXH LVTPGIQDSQ KVVIRAGISA GDRVVTDGID RLTEGAKVEV VEAQSATTPE EKATSREYAK KGARS
<b>Source</b>	Yeast
<b>Target Names</b>	mdtA
<b>Protein Names</b>	Recommended name: Multidrug resistance protein MdtA Alternative name(s): Multidrug transporter MdtA
<b>Expression Region</b>	22-415
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.