



Recombinant Escherichia coli Uncharacterized protein yfdF (yfdF)

Product Code	CSB-MP305173ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P76505
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain K12)
Purity	>85% (SDS-PAGE)
Sequence	MLPSISINNT SAAYPESINE NNNDEVNGLV QEFKNLFGNGK EGISTCIKHL LELIKNAIRV NDDPYRFNIN NSSVTYIDID SNDDTHITIG IDNQEPIELP ANYKDKELVR TIINDNIVEK THDINNKEMI FSALKEIYDG DPGFIFDKIS HKL RHTVTEF DESGKSEPTD LFTWYGKDKK GDSLAIVIKN KNGNDYLSLG YYDQDDYHIQ RGIRINGDSL TQYCSENARS ASAWFESSKA IMAESFATGS DHQVVNELNG ERLREPNDVF KRYGRAIRYD FQVDDAKYKC DHLKEIVSTL VGNKINVGHS QKIYKHFKDL EGKIEERLQN RQAEYQNEIN QPSAPGVNFD DI
Source	Mammalian cell
Target Names	yfdF
Protein Names	Recommended name: Uncharacterized protein yfdF
Expression Region	1-352
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.