



# Recombinant Escherichia coli Uncharacterized J domain-containing protein djIB (djIB)

<b>Product Code</b>	CSB-EP305198ENV
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P77381
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain K12)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MKNCWKILDI EETTDVDIIR RAYLALLPSF HPETDPQGFK QLRQAYEEAL RIAQSPAKSV WQPEEYEVAE HEILLAFRAL LASDSEFLP SAWQRFIQQL NYCSMEEIDE LRWSLCTIAM NTAHLSFECV VLLAERLRWL QEENTGEIDE EELESFLYAI AKGNVFNFT ILHLPVAVQN DTIDFYQMFA RIWSSHPQWL TLYLAQHRAV IIPDDAKLHR NLLRWYSAGR LDIPPELLDYA QSWRETEPDN EDAPYYEYAAQ RVYCGEGESL LAELCDYWRE YPSTQADALM LQWCRQHRVD YYPLLVMIE ARDLVNDQGK PLYVPGDSA RTRFHLYEIL SDEKLSALGR SLVEMVLHKG RKPRISLTRD TEHTLWPLYL VAKQLVQACQ PTEESLMPIV SRLDAENRCP LEALIIRLL IQAANFTEKQ TVEPEPQPQP MPVDDGGPGC LGIKIIFYI FIFAGLIGKI LHLFG
<b>Source</b>	E.coli
<b>Target Names</b>	djIB
<b>Protein Names</b>	Recommended name: Uncharacterized J domain-containing protein djIB
<b>Expression Region</b>	1-475
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.