



# Recombinant Methylaspartate mutase E chain (glmE)

<b>Product Code</b>	CSB-EP305262CLT-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P80077
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Clostridium cochlearium
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MELKNKKWTD EEFHKQREEV LQQWPTGKEV DLQEAVDYLK KIPAEKNFAE KLVLAKKKGI TMAQPRAGVA LLDEHIELLR YLQDEGGADF LPSTIDAYTR QNRVDECENG IKESEKAGRS LLNGFPGVNY GVKGCRKYLE AVNLPLQARH GTPDSRLLAE IIHAGGWTSN EGGGISYNVP YAKNVTIEKS LLDWQYCDRL VGFYEEQGVH INREPFGLT GTLVPPSMSN AVGITEALLA AEQGVKNITV GYGECGNMIQ DIAALRCLEE QTNEYLKAYG YNDVFVTTVF HQWMGGFPQD ESKAFGVIVT ATTIAALAGA TKVIVKTPHE AIGIPTKEAN AAGIKATKMA LNMLEGQRMP MSKELETEMA VIKAETKCIL DKMFELGKGD LAIGTVKAFE TGVM DIPFGP SKYNAGKMMP VRDNLGCVRY LEFGNVPFTE EIKNYNRERL QERAKFEGRD VSFQMVIDDI FAVGKGRLIG RPE
<b>Source</b>	E.coli
<b>Target Names</b>	glmE
<b>Protein Names</b>	Recommended name: Methylaspartate mutase E chain EC= 5.4.99.1 Alternative name(s): Glutamate mutase subunit epsilon
<b>Expression Region</b>	1-483
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.